HITACHI

U-3900/3900H



Hitachi Spectrophotometer U-3900/3900H Hitachi High-Technologies

Spectrophotometer Meeting a Wide Range of Analytical Needs from Liquid to Solid Sample Measurements

- Measurable over a broad absorbance range thanks to low stray light and low noise. (Model U-3900: -3.8 to 3.8 Abs, 0 to 300%T Model U-3900H: -5.5 to 5.5 Abs, 0 to 300%T)
- Stable monochromator in double beam optics (Baseline flatness Model U-3900 : within ±0.0003 Abs, Model U-3900H : within ±0.0004 Abs)
- Simple instrument control and diversified quantitative analysis supported by UV Solutions program for U-3900 (in connection with PC)
- A full range of accessories for covering both liquid and solid sample measurements



Stray light : 0.015% or less Photometric range : -3.8 to 3.8 Abs

> Two types available for selection according to sample and application purpose. Usable over an extensive field including analyses of water quality, the environment, biotechnology, drug manufacture and materials.



Stray light : 0.00025% or less Photometric range : -5.5 to 5.5 Abs

Adoption of Stigmatic Concave Diffraction Grating

Hitachi Model U-3900/U-3900H spectrophotometer adopts a Seya-Namioka mount monochromator and a stigmatic concave diffraction grating.

Because a concave diffraction grating is usable for both converging and dispersing light, it allows composition of an optical system with a small number of mirrors.

In this design, loss of light and aberration are suppressed, so a bright optical system can be configured.



Stable Optics with Double Beam



As a light source, a WI lamp (visible region) and a D2 lamp (ultraviolet region) are provided for selective use according to measuring wavelength range.

Double beam optics is adopted for ensuring stable measurements, in which the monochromatic beam selected with a monochromator is split into reference beam and sample beam with a rotating mirror (sector mirror) and the beams are directed into the sample compartment. In one model, the U-3900, a spherical mirror is used before the entrance slit. In the other model, the U-3900H, a grating is used before the entrance slit.

Since the Model 320 was launched in 1979, Hitachi medium-size spectrophotometers have been employed by customers in 25 countries.



Hardware

Hardware structure with priority given to ease of operation and data reliability.

Hardware structure attaching greater importance to ease of operation

USB communication is adopted between the spectrophotometer and PC. And, because the top face of the spectrophotometer is flat, a notebook PC can be mounted on it. Therefore, the spectrophotometer and PC can be connected promptly.





Incorporation of double monochromator

Due to mounting of a double monochromator which uses Hitachi's original stigmatic concave diffraction grating, an excellent linearity is ensured up to high concentrations. Hence, highly reliable quantitative analysis is possible.



Measurement with change in scan speed for ultraviolet region

Scan speed is changeable for the ultraviolet region. In this wavelength region, noise can be reduced by slowing down the scan speed. Owing to this feature, a noise-suppressed spectrum is obtainable over the entire range from visible to ultraviolet region by a single scan.



Effective in trace sample measurement

Satisfactory measurement is achievable even with 5, 25 and 50 μ L micro-sample cells because the beam is finely converged in the sample compartment.

Shown here are spectra determined in the ultraviolet region by measuring nucleic-acid adenosine with a micro-sample cell (internal volume $25 \ \mu$ L). A high S/N ratio was obtained.

Model : U-3900 Scan speed : 300 nm/min Slit : 2 nm



Ease of maintenance (in lamp replacement)

Lamp cable is connected by means of a socket, so each lamp can be removed or inserted without using a tool such as flat-head screwdriver.



Stable baseline

The Model U-3900 series assures a stable baseline in a wavelength range from 190 to 850 nm. (Baseline flatness Model U-3900 : within ± 0.0003 Abs, Model U-3900H : ± 0.0004 Abs)

Data can be measured stably even in a long-time measurement of enzyme activity, etc.





Model : U-3900 Slit : 2 nm Wavelength : 500 nm

Original differential feedback system

Sample signal, reference signal and zero point rise are always monitored and photomultiplier voltage is changed so that the sample or reference signal, whichever larger, becomes constant, whereby minus absorbance can be measured. Also, measurement in a broad dynamic range is allowed, e.g., difference spectrum measurement with different samples set on reference and sample sides.

> Model : U-3900H Scan speed : 300 nm/min Slit : 2 nm





Software

"UV Solutions for U-3900" program has been prepared for efficient instrument control and various quantitations.

A series of operations from analysis method setup to data processing can be initiated by clicking each button.





Analysis method setup window

Clicking this button enables the user to set analytical conditions such as measurement mode, measurement range and scan speed.

tem Series	 	
Sector State	 Net No. 1	
-		
The base late.		



Baseline measurement

Baseline measurement procedure can be started by clicking this button.

Upon measurement, data after baseline correction is obtainable.



Sample information setting

Clicking this button allows setting of detailed information about a sample to be measured, data saving location, etc.

Cowant	NV have served	L
VANTE	The second	
	1" Auto Text tim	
	(** Autor Marte Mile	
Fienane	Salert.	



Report Reporting

Data file can be output in a report format by clicking this button.

Enriched functions such as data comparison and preview are supported by UV Solutions for U-3900.

Easy comparison of measured data

Measured data can be compared easily by overlaying spectra or in Abs value at the specified wavelength. (A maximum of 10 spectra can be compared at 12 specified wavelengths).



Factor of data processing (quantitation) changeable

"Correction factor," "decimal digits of concentration" and "concentration unit" are settable on the sample table window. Setting can be determined in consideration of a sample to be measured, its concentration, etc.

New Caller	N Tarahaman Chara	Desen Talas
L	Hogen Linearthy Alarce Linearthy Linearthy 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000	
in the second se	jere i service service i grand interview i formation large tree i formation large tree i grand interview i formation grand interview i grand interview i grand interview i grand i grand interview i grand i grand interview i grand i grand i grand i grand i grand i grand i grand i grand i gra	ter ter p
l Suttine	8	



On the data processing window, factors such as "decimal digits of concentration" can be set.

10

File loading function with preview

Using a tool button for file loading with preview, max. two sets of data can be previewed without opening data files.

Previewable data loading button



Clicking a data file will display

the spectrum of the data (2 sets of data).



Reuse of analysis method for measured data

When it is desired to carry out measurement by the same analysis method as used for the already measured data, the "Apply analysis method" button is usable. The analysis method can be loaded and applied to a new measurement by clicking this button.



Control of lamp ON time

Total operation (ON) time of the WI and D2 lamps used in the U-3900 series can be checked on the software. This time counting is usable as a reference for judging the replacement time point for each lamp.

D2 Lamp	0	- h	Reset D2 Lamp
WE Lamp	0	h	Reset WI Lamp

Measured data exportable to commercially available software

Data such as measured spectrum can be pasted to Microsoft® Word and Microsoft® Excel, and converted into an ASCII text file. Using such software, a report form can be edited.

Page Report	t data laite		a)
		Include data living	
 Provide mellical Provide mellical Provide graph Provide data bring Provide parts total 	P Start VL Control D Value VL Control D Value VL Value Start	Per Alber	-
Face Face - Face range - MS Serv S	ed .		
Perfilme 30	- peec	18 E	

Accessories Expediting Application to Multi-Sample Measurement, Micro-volume Sample and Many Others

Micro cell holder

(P/N 122-0060)

Suitable for micro-sample measurement in medical and biochemical fields.

Specifications			
Wavelength range	220 to 800 nm		
Repeatability in cell placement	Within ±0.3%T		
Baseline flatness	Within ±0.0025 Abs (when using 50 µL micro-sample cell)		

Micro-sample cell

The following cells are usable for the above-mentioned micro-sample cell holder.

Part No.	Part name	Capacity	Optical path length
130-0622	50 µL cell	50 µL	10mm
130-0623	25 µL cell	25 µL	5mm
130-0621	5 µL cell	5 µL	0.5mm

Auto sipper (P/N 2J1-0105)

This computer-controlled sample sipper is provided with a sample recovery function and other versatile functions. In combination with an autosampler, this unit makes it possible to carry out automated labor-saving analysis.

Specifications			
Minimum sample volume	0.6 mL		
Carryover	1% or less		
Gell capacity	Approx 60 µL		
Ciptical path length	.10 mm		
Reference beam side	10 mm rectangular cell mountable		

Electronic thermostatted auto sipper (P/N 2J1-0106)

The flow cell section is maintained at a constant temperature level under accurate control.



Specifications (reference beam side not temperature-controlled)

Minimum sample volume	0.6 mL
Carryover	1% or less
Cell capacity	Approx. 50 µL
Optical path length	10 mm
Setting temperature	20 to 40°C
Setting accuracy.	Within ±0.5°C
Reference beam side	10 mm rectangular cell mountable

AS-1010 autosampler (P/N 2J1-0121/0122)

This unit is used for multiple-sample measurement in combination with an auto sipper or in flow injection analysis. A suction needle can be moved in three directions X, Y and Z.



Specifications (sample tube not included)
Sample tube size
Outside dismeter 15 mm, height 105 mm (option required)
Outside dismeter 12 mm, height 105 mm

Water circulated cell holder (P/N 210-2111)

Water from a thermostatic oven is circulated through this cell holder to maintain a sample cell at a constant temperature. (Temperature control : R and S)



 Specifications (circulatory thermostatic oven and cell not included)
 Operating temperature range

 Operating temperature stability
 Room temperature to 40°C

 Temperature stability
 Within ±0.3°C

Electronic thermostatted cell holder (P/N 131-0306/0307)

This cell holder comes standard with an incorporated magnet stirrer. The temperature of sample in a cell is maintained at a constant level, and a temperature value down to 0.1°C can be indicated. Since this unit is of electronic thermostatted type with forced air cooling, quick heating and cooling can be performed without a water circulating thermostatic oven. (Temperature control : S only)



Specifications

inder condition of 25°C room temperature)
Vithin ±2°C (*) (difference between set emperature and actual sample temperature).
Vithin ±0.5°C (*)
0 mm dell (dell not supplied)

Electronic thermostatted cell holder (P/N 131-0301/0302)

In protein and nucleic acid melting measurement, sample temperature can be changed continuously to determine variation in absorbance. Being of an electronic thermostatted type, this cell holder is capable of quick heating and cooling. Sample temperature can be increased and decreased isothermally. Because this holder is equipped with a stirrer, the internal cell temperature can be kept uniform. (Temperature control : R and S)



 Specifications (circulatory thermostatic oven not included)

 Applicable cell
 10 mm cell (not included in this unit)

 Temperature range
 0° to 100°C (settable in increments of 0.1°C)

 Temperature control accuracy
 Within ±2°C (*) (difference between set temperature and actual sample temperature)

 Temperature stability
 Within ±0.5°C (*)

 Provided with an isothermal regulating function

 Room temperature : 25°C, sample : distilled water, circulatory water temperature : 22°C

Setting temperature : 10°C to 60°C

A circulatory thermostatic oven needs to be prepared separately.

Micro flow cell unit (P/N 210-2113)

Suitable for continuous measurement of a micro-quantity of sample.



Specifications	
Cell capacity	70 µL
Optical path length	10 mm (quartz flow celi used)
Connection tubing	Tetlori tube of outside diameter 2 mm and inside diameter 1 mm

Flow cell unit

(P/N 210-2173)

The inside of this cell is structured to minimize stagnation of liquid and adhesion of air bubbles.



Specifications

Cell capacity	600 µL
Optical path length	5 mm (quartz flow cell used)
Connection lubing	Teflon tube of outside diameter 4 mm and inside diameter 3 mm
Reference beam side	5 mm rectangular cell (standard accessory)

LC flow cell unit

(P/N 210-2131)

A flow cell especially designed for liquid chromatography.



Specifications		
Cell capacity	8 µL	
Optical path length	8 mm (quartz flow cell used)	
Baseline flatness	±0.001 Abs (200 to 350 nm) +0.002 Abs (190 to 850 nm)	

6-cell positioner with temperature control (P/N 2J1-0103/0104)

Six 10 mm cells can be mounted on the sample beam side, and they can be changed over automatically at certain intervals. (Temperature control : S only)



Specifications			
Repeatability in cell changeover	Within ±0.5% (at 100%T)		
Applicable cell	10 mm cell (not included in this unit)		
Setting temperature	20 to 40°C		

Tandem cell holder

(P/N 210-2115)

A maximum of three 10 mm cells can be mounted on each of the sample and reference beam sides. Sample temperature can be maintained at a constant level by circulating temperature-regulated water through the cell holder section. (Temperature control : R and S)



Specifications (not including circulatory thermostatic oven and cell)
Temperature range
15 to 40°C
Temperature stability
+0.3°C

4-position rectangular long-path cell holder

(P/N 150-0940)

Four rectangular long-path absorption cells can be mounted on the sample beam side, and they can be changed over externally.

Specifications

Cell length 100 mm, 50 mm to 10 mm cells applicable

5-position turret cell holder (P/N 210-2110)

Five 10 mm rectangular cells can be mounted on the sample beam side, and a micro-cell mask (200-1537, 200-1538) can be inserted in each cell holder. (Cells and micro-cell mask are not included.) It is recommended to prepare a set of five cells.



Rectangular long-path cell holder (P/N 210-2107)

Rectangular cells having the following optical path lengths are applicable: 10, 20, 30, 40, 50 and 100 mm. The outer width is 12.75 mm.



Rectangular cell

(Note) Absorption cells are not included in the main unit.

an No.	Fan name	
123-1004	10 mm quartz cell set (2 cells in set)	
123-1010	10 mm glass cell set (2 cells in set)	
124-0349	30 mm quartz cell set (2 cells in set)	
124-0350	50 mm quartz cell set (2 cells in set)	
210-3939	100 mm guartz cell set (2 cells in set)	

Examined 10 mm rectangular quartz cell

Optical path lengths measured at 21 points on cell using a three dimensional measuring instrument are indicated down to the fourth decimal place (mm).

Part No. Part name 210-1462 Examined 10 mm quartz cell (1 cell)

Mask for micro cell

For use of this mask, it is inserted in a 10 mm rectangular cell holder which is standard with the main unit.

Part No. Remarks

Cell

200-1537 Mask width 1.5 mm 200-1538 Mask width 1.2 mm

 Part No.
 Part name

 124-0357
 Micro 10 mm quartz ceil (2 ceils in set)

 200-0551
 Black micro 10 mm quartz ceil (2 ceils in set)

Cylindrical long path cell holder (P/N 210-2108)

Applicable cylindrical cells are listed below.



Applicable cylindrical cells

23-1023	10 mm quartz cylindrical cell for liquid sample (2 sells in set)
23-1024	10 mm glass cylindrical cell for liquid sample (2 sells in set)
23-1027	20 mm quartz cylindrical cell for liquid sample (2 sells in set)
23-1028	20 mm glass cylindrical cell for liquid sample (2 sells in set)
23-1031	50 mm quartz cylindrical cell for liquid sample (2 sells in set)
23-1032	50 mm glass cylindrical cell for liquid sample (2 sells in set)
23-1035	100 mm quartz cylindrical cell for liquid sample (2 sells in set)
23-1036	100 mm glass cylindrical cell for liquid sample (2 seils in set)
23-1043	10 mm quartz cyfiridrical cell for gas sample (2 sells in set)
23-1044	10 mm glass cylindrical cell for gas sample (2 sells in set)
23-1047	20 mm quartz cylindrical cell for gas sample (2 sells in set)
23-1048	20 mm glass cylindrical cell for gas sample (2 sells in set)
23-1051	50 mm quaitz cylindrical cell for gas sample (2 sells in set)
23-1052	50 mm glass cylindrical cell for gas sample (2 sells in sel)
23-1058	100 mm quartz cylindrical cell for gas sample (2 sells in set)
23-1056	100 mm glass cylindrical cell for gas sample (2 sells in sel)

Glass filter holder

(P/N 210-2109)

Used for transmittance/absorbance measurement of a solid sheet sample such as glass filter.



Specifications

Sample Inickness	0.5 10 5 mm	
Sample size	Minimum : 12 x 25 mm Maximum : 55 x 100 mm	

Film holder



Convenient for measurement of film-shaped samples.



φ150 integrating sphere accessory (P/N 2J2-0175)

Designed for diffuse reflectance measurement of a solid sample surface and absorbance measurement of a turbid sample. With an aperture ratio as small as 2%, this unit is usable for high-accuracy colorimetric measurement.



Specification

350 to 750 nm	
±0.5%T	
2%	
Mountable	
	350 to 750 nm ±0.5%T 2% Mountable

φ60 integrating sphere accessory (P/N 2J2-0176)

Designed for absorbance measurement of a turbid sample and reflection measurement of a solid sample surface.



pecifications

opeonition	
Wavelength range	250 to 800 nm
100%T line flatness	±1%T
Aperture ratio	7.8%
Specular reflection measurement attachment	Standard-equipped

5° specular reflectance accessory (P/N 2J2-0177)

Using mirror reflection of a sample, relative reflectance is measured with respect to the standard reflection plate (aluminum-evaporated plane mirror). Applicable to film thickness measurement and spectral reflectance measurement.



Specifications Angle of incidence 6* Sample area 25 mm in diameter or more

Polarizer holder (P/N 210-2130)

Sample beam is linearly polarized for measurement of polarization characteristics or a sample is placed between the polarizer and analyzer for measurement of optical rotary power.



Specifications		
Wavelength range 400 to 750 nm		
Sample area	Minimum 12 mm x 25 mm Maximum 55 mm x 100 mm	
Sample thickness	0.5 to 5 mm	

Optional program

-		
Part No	Part name	
2.12-0317	GLP/GMP program	

U-3900/3900H

Model	U-3900	U-3900H	
Monochromator	Diffraction grating Single monochromator Seya-Namioka mount	Diffraction grating-diffraction grating Double monochromator Seya-Namioka mount	Spectrophotomet
Wavelength range	190 to 9	00 nm ^(*1)	control
Spectral bandpass	0.1, 0.5, 1, 2, 4	4, 5 nm (6 steps)	
Stray light	0.015% (Nal : 220 nm. 1	0.00025%	
Wavelength accuracy	±0.1 nm (at 656.1 nm af	iter wavelength calibration)	
Wavelength setting repeatability	±0.0	05 nm	M
Photometric mode	Abs, %T, %	R, E(S), E(R)	conditions
Photometric range	Abs : -3.8 to 3.8 Abs (effective range) %T : 0 to 300%T	Abs : -5.5 to 5.5 Abs ^(*2) (effective range) %T : 0 to 300%T	Conditions
Photometric accuracy (checked with NIST SRM930)	±0.002 Abs (0 to 0.5 Abs) ±0.003 Abs (0.5 to 1.0 Abs) ±0.006 Abs (1.0 to 2.0 Abs) ±0.3%T		Execution of
Photometric repeatability (checked with NIST SRM930)	±0.001 Abs (0 to 0.5 Abs) ±0.0015 Abs (0.5 to 1.0 Abs) ±0.003 Abs (1.0 to 2.0 Abs) ±0.1%T		measurement
Response	High resolution, Standard		
Baseline flatness	Within ±0.0003 Abs (190 to 850 nm)	Within ±0.0004 Abs (190 to 850 nm)	
Baseline stability	Within 0.0002 Abs/hr (at 50	0 nm, 2 hours after power-on)	Departing/
Baseline memory	3 channels (system : 1 c	hannel, user : 2 channels)	Display
Wavelength scan speed	1.5, 3, 15, 30, 60, 120, 300, 600, 1200, 1800, 2400 nm/min		Diopidy
Light source	Adjustment-free deuterium lamp D2 lamp): Ultraviolet region Adjustment-free tungsten iodine lamp (50 W)(WI lamp) : Visible region		
Light source changeover	Automatic changeover linked with wavelength Changeover wavelength : Selectable in a range of 325 to 370 nm)		
Sample compartment	Beam spacing: 100 mm 120 (W) × 300 (D) × 140 (H) mm		
Detector	Photon	nultiplier	
Data processing unit	PC: OS Windows [®] XP Professional SP2		
Dimensions (spectrophotometer main unit)	680 (W) × 692 (D) × 257 (H) mm		Data processing
Operating temperature/humidity	Temperature : 15 to 35°C. Humidity: 25 to 85% (non-condensing)		
Weight	45 kg (spectrophotometer main unit)		
Power consumption	100 V AC 50/60 Hz, 300 VA (excluding personal computer and printer)		
UV Solutions program	Standard	software	

Software Functions (common to U-3900/U-3900H)

	Wavelength/Time Scan, Measurement and Data Processing		Photometry
Spectrophotometer control	 Wavelength shift (Go to λ) 100% T adjustment (auto zero) Automatic wavelength calibration Detector zero adjustment 		
Moosuring	Measuring condition setting Condition loading Condition saving (desired number of files, file overwriting/deletion possible) Automatic start function (measuring conditions automatically set upon startup of software)		
conditions	_		Condition setting for working curve (1st to 3rd order, segmented line) Standard data setting (20 standards, average of 20 data values) K factor input
Execution of measurement	 Measurement of spectrum/change w Repetitive spectrum measurement Setting of sampling interval Measurement with scan speed c in ultraviolet region 	ith time hanged	Remeasurement of working curve Interrupt measurement
	 Baseline measurement (3 chann (system baseline : 1 channel, us) 	els) er baselir	ne: 2 channels)
	 Sample name Comment input Ruled line printout ON/OFF 	ut Measurir	ng condition printout ON/OFF
Recording/ Display	 Printout and display of spectrum change with time Spectrum loading (with preview) Spectrum saving Spectrum printout/display 	/	Working curve printout/display Data deletion Data loading Data saving Data list printout/display
Data processing	Rescaling (numerical value input, cursor input) Spectrum trace Smoothing Peak detection Data printout (printout of wavelength/time at fixed intervals) Graph axis conversion Spectral calculation (arithmetic calculation/coefficient calculation) Differentiation (1st to 4th order) Data reset Rate calculation (only in time scan mode) Spectrum display window		Working curve trace Working curve erasure Data printout Sample data erasure Statistic calculation Calculation of determination coefficient Setting of correction coefficient Setting of number of decimal places for concentration Setting of concentration unit
Miscellaneous	Viscellaneous File conversion (ASCII/JCAMP) Setting of number of decimal places for display Cell length conversion Data transfer/graph copy to Microsoft [®] Excel Print preview Display of lamp ON time		
S	Standard Equipment Q'ty		

 Standard Equipment
 Q ty

 Spectrophotometer main unit
 1 set

 Tools
 1 set

 Instruction manual
 1 set

NOTES: 1. Absorption cells are not included in the standard equipment, and thus should be prepared separately. 2. A PC set is not supplied as standard equipment. It should be prepared separately.

Notice: For proper operation, follow the instruction manual when using the instrument.

Specifications in this catalog are subject to change with or without notice, as Hitachi High-Technologies Corporation continues to develop the latest technologies and products for our customers.

@Hitachi High-Technologies Corporation

Microsoft® Excel, Word and Windows® XP are trade names of Microsoft Corporation in U.S.A.
 Microsoft® and Windows® are registered trademarks of Microsoft Corporation in U.S.A.

Tokyo, Japan

Within ±0.0004 Abs (U-3900H) *2: With 1%T light attenuating plate

http://www.hitachi-hitec.com/global/science/

24-14 Nishi-Shimbashi 1-chome, Minato-ku, Tokyo, 105-8717, Japan Tel: +81-3-3504-7111 Fax: +81-3-3504-7123

PRINTED WITH SOY INK For further information, please contact your nearest sales representative.